

BEAUCCLARKE

For Private Circulation Only.

Chas. Beaucclarke

(PROP. J. BLOCK).

**6, NEVILLE STREET,
(opposite Central Station)
NEWCASTLE-ON-TYNE.**

**BIRTH CONTROL
SPECIALISTS.**

INTRODUCTION.

To most people the question of limitation of family is, from a variety of causes, of the highest importance, and to many a matter of absolute necessity, the only question being as to the safest and most convenient method to adopt. In this handbook you will find various methods of Birth Control. With regard to choice of method this must be left to the individual requirements of each couple.

It is most important that all Birth Control Appliances should be of the most reliable hygienic make and only proved brands used. We have been established since 1894, and our goods are known throughout the United Kingdom for their quality and reliability.

RUBBER SHEATHS **(Method No. 1.)**

The Sheath is the most commonly used method of Birth Control and is a safe and sure check.

The Sheath is composed of fine rubber closed at one end and is worn by the male. They are simply rolled on. It is advisable to leave the closed end a little loose. Teat-end sheaths have a receptacle on the closed end which minimises the danger of breakage. Do not wash ordinary sheaths. If you wish to wash them, then purchase the washable brand.


X
✓

Beaclarke's have been Established since 1894.

Every sheath supplied by us is subjected to the most severe test before leaving the works, so that the slightest flaw may be detected. It is, therefore, evident that no sheath with the least fault could pass such a test (applied to every quality of sheath supplied by us.) The dearer grades of sheaths are thinner and almost undetectable, but entirely without possibility of leakage.

All our Sheaths bear the famous "**Dreadnought**" Brand and are only obtainable at our **Newcastle** address; beware of imitations.

We have no connection with any other firm trading under the same name.



Chas. Beaclarke
PROP. J. BLOCK.

G. NEVILLE STREET,
NEWCASTLE-ON-TYNE.

STRONGEST, CHEAPEST & BEST.

All Sheaths are rolled and packed in envelopes with the above design printed on—refuse any without.

Beaclarke's for Hygienic Birth Control Appliances.

Beaclarke's Famous "Dreadnoughts."

Seamless Silk-finished Preventatives

6d. each, **2/6** half doz., **5/-** doz.

Finest Ivory Seamless Silk-finished

Preventatives, 1st quality.

1/- each, **5/-** half doz., **10/-** doz.

Finest Ivory Seamless Extra Stout Silk-

finished Preventatives, Finest quality,

1/6 each, **7/6** half doz., **15/-** doz.

Special "Transparent Dreadnought"

The strongest and thinnest Preventative made. It has taken many years to produce this Sheath to perfection, but now it is the article *par-excellence* of exceptionally thin, fine silky texture and elasticity.

1/6 each, **7/6** half doz., **15/-** doz.

BOX OF THREE—An Aluminium Case containing 3 seamless silk-finished Preventatives **2/-**.

The Paragon Washable Sheath.

These are a stout-made Preventative of the finest quality and prepared to last a considerable time thus proving a great saving to persons of limited means.

2/-, 2/6 and 3/6.

Beauclarke's for Proved Brands.

The Poor Man's Friend.

This Preventative is manufactured from pure medicated cut sheet rubber and prepared in such a way that it can be used any number of times. Will last for years. This Sheath can also be rolled up and used as a Check Pessary.

5/6 each.

We also manufacture a special Washable Sheath, as above, with a moulded double top, which we can specially recommend.

6/6 each.

"Beauclarke's" Seamless Transparent Sheath

This is the lightest and strongest washable sheath sold to day. It has taken many years to get this sheath perfect, and be able to sell at the extremely low figure of

4/6 each.

This is our own **Special Line** in Washables.

"Monarch" Seal Brand (Regd.)

These are a transparent seamless sheath, lighter than the ordinary Paragon sheaths, but very strong and will last for years. We can highly recommend this sheath.

7/6 each.

Spanish Skins (also known as Fish Skins).

These are made of a very thin membrane, fine and close fitting and contain no rubber whatever. After fitting, damp with water and they will conform naturally to shape.

1/6 each, 7/6 half doz., 15/- doz.

Beauclarke's known throughout the United Kingdom.

American Tips (or Short Sheaths).

American Sheaths, also known as Malthus Caps, are short sheaths for use where complete covering is not desired and are fitted with a strong elastic band which prevents the slightest possibility of them coming off.

8d. each, 3/6 half doz, 7/- doz.

Soluble Quinine Pessaries (Method 2)

The Wife's Friend,

The most successful and simple method of introducing into the Vagina that powerful Contraceptive Quinine and Cocoa Butter. They are made in the form of small flat dome-shaped cones which can be conveniently and unobtrusively slipped into the Vagina a few moments before intercourse.

They dissolve and set free an ingredient which instantly destroys the vital properties. These are all treble strength and absolutely the strongest obtainable.

2/6 per box (12 in box).

"MONARCH" BRAND (Regd.)

Extra strong perfumed Pessaries **3/6 box (12 in box).**

CHINOSOL PESSARIES.

These Pessaries have the same action as Quinine Pessaries and are recommended by Marie Stopes in "Wise Parenthood."

2/6 box of 12.

"Q" CONES (Red Seal) Quinine Pessaries.

These are a non greasy and odourless Pessary and are now being used in place of the greasy Quinine Pessary. Absolutely reliable.

3/- per box (12 in box),

Beaclarke's—The Reliable House.

Anti-Conception Check Pessary (Method 3.)

A simply devised appliance of pure medicated rubber. It is constructed on a commonsense principle, can be worn any length of time with ease and comfort, is easily adjusted and removed. This Preventative is worn by the woman. It is dome-shaped in construction and when fitted covers the neck and mouth of the womb during intercourse. It is advisable to use a Quinine or Chinosol Pessary with the Check Pessary.

They are made in three sizes. No 1 for small-built women without children or one child. No. 2 for average built women with 2 or 3 children. No. 3 for large built women with 3 or more children.

Solid Rim, best quality 3/6 each.

„ 2nd „ 3/ „

Contraception Sponge (Method No. 4.)

This is a small round prepared sponge incased in a silk net with cord attached. They are moistened with water, rubbed with compound quinine ointment or dipped in Compound Quinine Solution and inserted into the Vagina as far as possible. As a double precaution a Soluble Quinine of Chinosol Pessary should also be inserted. It is safest to leave the sponge in until the next morning.

Prices 1/- each, 10/- doz., (Medium size.)

„ 1/3 „ 12/- „ (Large size.)

Douches (Whirling Spray). Method No. 5.

This method consists of syringing the Vagina immediately after intercourse. The solution for the syringe is made by dissolving a Compound Quinine Powder in water.

Beaclarke's give Advice FREE.

Douches—continued.

Every woman should use a Whirling Spray periodically as it is well known that regular injections are recommended very strongly by doctors. They recognise its great value, and most of the common ailments to which the great majority of women are subject could be easily prevented by regular use of a

WHIRLING SPRAY.

This Spray is also very useful as a further additional safeguard after the use of Quinine, Chinosol or Check Pessaries, also after using a safety sponge.

No woman should be without one. Made of the finest British-made rubber and vulcanite with extra large bulb.

6/6, 7/6, and 8/6.

“SUPER” WHIRLING SPRAY.

This is a very high-class spray with a powerful spraying action and extra large bulb. 10/6 each.

SPRAY TABLETS.

These are specially compounded for use with Whirling Spray or Enema.

Box of 36 Tablets 1/9.

QUININE COMPOUND POWDER.

These are also for douching purposes and are a reliable preventative of conception.

2/- per doz. Powders.

Advice Free and without Obligation.

ENEMA SYRINGE.

We have a very large selection of these Syringes from
2/6, 3/6, 4/6 and 5/6.

SUSPENSORY BANDAGES.

2/6, 3/6, 4/- and 5/6 small, medium or large.

JOCK STRAPS.

4/6 each, small, medium or large.



We specialize in
**SURGICAL BELTS, TRUSSES. ELASTIC
HOSIERY, ARTIFICIAL LIMBS,
CRUTCHES, BED PANS, URINALS,
RUBBER URINALS** for day and night use,
and every description of Surgical Appliances.



MARIE STOPE'S EDITIONS

OF

"Married Love," 6/-. "Wise Parenthood," 3/6.
"Radiant Motherhood," 6/- and "Contraception," 15/-.

KEPT IN STOCK.

Instructions to Purchasers.

We respectfully request customers in ordering to give full name and address and to write clearly. State the number of page on which article appears.

When ordering Sheaths state whether **Plain** or **Teat** are required.

All orders should be accompanied by Postal Order, Money Order or Cheque, and should be crossed thus & Co. and made payable to Chas. Beauclarke.

All goods are forwarded under plain cover.

Orders of 2/6 and over, Carriage paid, excepting Marie Stope's Publications. These are all 6d. extra for postage above catalogue prices.

BEAUCCLARKE

FOR PRIVATE CIRCULATION ONLY.

PRICE LIST

OF HIGH-CLASS

Surgical Rubber, Medical
and Surgical Appliances.

NOTE OUR ONLY ADDRESS:—

CHAS. BEAUCCLARKE

(Proprietor—J. BLOCK),

6, NEVILLE STREET,
NEWCASTLE-ON-TYNE.

FACING TYNEMOUTH ENTRANCE CENTRAL STATION.

SEND YOUR ORDER TO US.

ALL PARCELS FORWARDED UNDER PLAIN COVER.

ALL OUR GOODS ARE GUARANTEED.



Ladies Please Write c/o Manageress.