


PATENTEX D.R.P.

The ideal conception preventative.



What are the essentials in a conception preventative that it may be used with the utmost confidence?

1. It must be beyond question absolutely reliable und effective.
2. It must be non-injurious to health.
3. Its operation must be simple and clean.

Do the so called preventatives so largely advertised fulfil these conditions?

No! They do not.

The well known American and French specialities are strikingly characterised by the saying:
„They are armours against pleasure and only cobwebs against conception.”

None of them could be termed abso lutely safe, to say nothing of their being troublesome and u
esthetic in their use, and a danger to health and life.

The ever progressing science of our time has evolved a change. Years of study and untiring dili-
gence, based on newly discovered scientific facts, have united ALL necessary qualities in a most ideal
anner, in the product which is known under the name of „Patentex” (protected by Patent in all countries).

The astonishing simplicity of „Patentex” can be seen from the above illustration. The pipe has
only to be inserted into the female vagina and the thumbscrew given a quarter to half a turn; by so
loing the preventative substance will pass through the tube and reach its intended place in the mouth of
he womb.

The substance, a thick fluid, closes firstly the mouth of the womb so that the seed cannot reach
he female ovary; secondly it dissolves the albumen of the seed into shreds, and thirdly it kills the seed
germ ad the moment of effusion. The threefold protection, examined by experts to the smallest detail
and pronounced effective, makes conception impossible.

The preparation also possesses other qualities: it is clean, hygienic, beneficial. and health producing
to the female organs.

„Patentex“ possesses every desired quality. That is, besides being a certain and immediate preventative, simple and without danger in its use, it is beneficial to health, clean and cheap, and may be our fullest confidence. It offers no disadvantages whatsoever, and is without doubt the most IDEAL PREVENTATIVE AGAINST CONCEPTION, and the surest safeguard against the danger of infection.

Instructions for the use of „Patentex“

Put the different parts of the apparatus together as shown in the illustration, and turn the handle of the tube a little to ensure that the substance leaves the pipe. If this is the case, then insert the tube into the female vagina so far, that besides the tin tube, only a small piece of the Glass Pipe is visible. Now give the key very slowly a quarter to half a turn so that the preparation enters the vagina. Turn out the apparatus and the injected substance acts for the next few hours as a sure safeguard against conception.

Guarantee!

The manufacture of „Patentex“ takes place under constant supervision of an Apothecary and certified Chemist.

Further advantages:

1. Against Nervousness.

The nervousness and the quick decay of youth and beauty of many women are known to be a result from the continual fear of pregnancy. These nervous women can have, by using „Patentex“, the conviction of a complete surety and with this a guarantee for a restoration to health.

2. Syringing.


Instead of the usual syringing, which every woman who is considerate of her health performs, and which occupies a long time, one injects a little Patentex into the female vagina, leaves it here for an hour or better still, from night till morning and then rinses with a syringe and lukewarm water. This method is most hygienic and most effective.

„Patentex“ is protected in all civilized states by Patent, Letter Patent etc.

|| Sole rights reserved. Infringements will be prosecuted. ||

Directions for the use of Patentex.

Patentex replaces injections as a prophylaxis measure against Infection, Leucorrhoea, catarrh of the vagina, and all other affections.



Screw the rubberpiece on to the tube such as shown in above sketch and turn the key a little so to make sure that the substance is forced out of the opening. Insert the frontpiece into the vagina until only a small portion of it is visible. Turn the key slowly a quarter or half way around, so that the substance is forced into the interior, then remove the apparatus. Once screwed together the apparatus need not be taken apart again.

Patentex against Leucorrhoea, catarrh of the vagina etc.

In case your physician does not direct otherwise, use Patentex 2 or 3 times per week principally in the evening, followed by a luke warm water syringing the next morning. To effect a permanent cure it is only necessary to repeat the application from time to time.

Patentex as a prophylactic measure against contagious diseases.

One application of Patentex gives absolute guarantee against contagion for several hours.

Patentex for Syringing.

Instead of the ordinary syringing, which every woman, who values her health, should carry out and which should last for some time, insert a small portion of Patentex into the upper vagina and allow same to remain there for an hour, or still better, until the next morning, then give a syringing with warm water. This method is the most hygienic.

Important

In order not to interfere with coception, do not use Patentex immediately before coitus. On account of its anticonceptional qualities Patentex is used by physicians in cases, where a reliable Anticoncipiant is indicated.

Patentex is protected by patent rights in all civilized countries.

The preparation is made up under constant supervision of a government chemist and analyst, which guarantees absolute security and harmlessness in each individual case.